

Biotin-Interferenz führt zu falsch ...

... niedrigeren Messwerten:	
AFP	Anti-HAV IgM
CA 125	HBs-Ag
CA 15-3	HBe-Ag
CA19-9	Anti-HBs
CA72-4	Anti-HBc IgM
CEA	Anti-HCV
C-Peptid	HIV-Suchtest
Ferritin	Lues-Suchtest
Insulin	FSH
Myoglobin	HCG
NT pro BNP	LH
PCT	Prolactin
S 100	PSA, total
SCC	PSA, Frei
Troponin T	PTH
Holotranscobalamin	SHBG
	TSH

... höheren Messwerten:	
Tg-Ak (TAK)	Anti-HAV
Anti-TPO	Anti-HBc
TSH-Rez. AK (TRAK)	Anti-HBe
Cortisol	
DHEA-S	
Estradiol	
Folsäure	
FT3	
FT4	
Progesteron	
Testosteron	
Vitamin B12	
Vitamin D 25-OH	

Quelle: ECLIA Testanleitungen, COBAS, Roche
Stand: 29.07.19